




Images for illustrational purposes only
Sensor should be placed underneath the bed sheet

Bed/Chair Sensor Paper Thin™

Available exclusively to Local Authorities from Possum



The Possum Paper Thin™ Bed Sensor

Unlike most commonly used sensors, the Paper Thin™ mat detects an individual's dielectric signals NOT their weight.

Bed/Chair Occupancy and Falls Prevention

- Proven to significantly reduce false alarms
- Can be cut to size and positioned to provide the best level of protection
- Works with any type of bed/mattress, including profiling beds
- Easy to fit - place underneath the bed sheet and attach to the mattress with the strong adhesive tape (which is positioned on the reverse of the mat)
- Cheaper than your current solution? (Challenge us.....)
- Can be used with the majority of suppliers Dispersed Alarm Units (DAU), care phones, paging solutions
- 100% Waterproof
- Easy to clean
- Tear resistant
- Tissue viability approved and permeable, assisting in the prevention and management of pressure damage

Fall Prevention

Evidence shows that effective fall prevention as part of a multi-factorial fall prevention programme (MFFP) is cost effective, significantly reduces fall rates (biggest cause of death in over 65s in hospital following hip/knee replacement in the UK) and simple to install and operate.*

Available exclusively to the Telecare Market from Possum.

Product Code: 035-082

Replacement Bed Sensors (pack of ten) Product Code: 035-103

Possum Ltd. 8 Farmbrough Close Stocklake Park Industrial Estate Aylesbury, Buckinghamshire
HP20 1DQ, United Kingdom

T: +44 (0)1296 461000 F: +44 (0)1296 461001 E: telecare@possum.co.uk
W: www.possumtelecare.co.uk

Disclaimer: Images for illustrational purposes only and the actual product supplied may vary in appearance. Product functionality, features and quality will not be affected.

*Nursing Times 24.05.11 / Vol 107 Online Edition / www.nursingtimes.net